



Fluorescent Proteins and Antibodies


For research use only

Fluorescent Proteins

- CoralHue[®] Azami-Green
- CoralHue[®] Dronpa-Green
- CoralHue[®] Kaede
- CoralHue[®] Keima-Red
- CoralHue[®] Kikume Green-Red
- CoralHue[®] Kusabira-Orange
- CoralHue[®] Midoriishi-Cyan


Fluorescent Protein Antibodies

CoralHue[®] Azami-Green




Bright Green!

CoralHue[®] Keima-Red



Large Stokes Shift!

CoralHue[®] Kaede



Photoconverting! *

CoralHue[®] Dronpa-Green






Photo-Activation! **

CoralHue[®] Kikume Green-Red




Photoconverting!

CoralHue[®] Kusabira-Orange



Bright Orange! **

CoralHue[®] Midoriishi-Cyan




Bright Cyan!

* Photo provided courtesy of Dr. Miyawaki, RIKEN Institute, Japan

**Photos provided courtesy of Dr. Michael Davidson, National High Magnetic Field Laboratory, The University of Florida.

CoralHue[®] Azami-Green




CoralHue[®] Azami-Green is derived from the stony coral whose Japanese name is "Azami-Sango". The CoralHue[®] Azami-Green (AG) fluorescent protein absorbs light maximally at 492 nm and emits green light at 505 nm. AG is stable in both acidic and basic conditions and does not show a significant loss of signal, giving it an advantage over other fluorescent proteins such as GFP. AG also matures rapidly to form tetramers that are highly fluorescent. This allows AG to be used to identify cells or to report gene expression without problems stemming from protein aggregation. AG has also been engineered as a monomeric fluorescent protein which allows AG to be used in protein fusion and subcellular structure studies. AG's tight tetramers and monomeric form give AG an advantage over GFP, whose oligomeric form limits GFP to cell labeling.



CoralHue[®] mAG1 expression in HeLa cells


AG has been engineered to provide several humanized, monomeric forms of the Azami-Green fluorescent protein which are useful in cases where monomers are preferred over tetramers. CoralHue[®] Azami-Green is available in several different plasmids, including expression plasmids, which can help to customize your research. AG is available as a cDNA plasmid which can be manually inserted into a plasmid in order to tag particular proteins of interest. AG is also available in plasmids suitable to construct C-terminal or N-terminal fusion proteins. Finally AG is available in several targeted expression plasmids to label the endoplasmic reticulum, the nucleoplasm, the mitochondria, or the plasma membrane.

| | Excit. /Emiss.Maxima (nm) | Extinction Coefficient(M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|-------|---------------------------|--|----------------------------|-----------------------|
| mAG 1 | 492 / 505 | 55,500 (492 nm) | 0.74 | pK _a = 5.8 |



CoralHue[®] mAG1 Mitochondria

Endoplasmic reticulum



Plasma Membrane


Nucleoplasm

References

- Shaner, N.C., et al., (2005) Nat. Methods. 2, 905-909. Review.
- Karasawa, S., et al., (2003) J. Biol. Chem. 278, 34167-34171.

For more information, go to www.mblintl.com

Vector



| CoralHue® Azami-Green Fluorescent Proteins | | |
|--|---|-------|
| Code No. | Product | Size |
| AM-V0021 | CoralHue® Azami-Green (AG-S1) | 20 µg |
| AM-V0031 | CoralHue® Monomeric Azami-Green (mAG1-S1) | 20 µg |
| AM-V0032 | CoralHue® Monomeric Azami-Green (pmAG1-MC1) | 20 µg |
| AM-V0033 | CoralHue® Monomeric Azami-Green (pmAG1-MN1) | 20 µg |
| AM-V0034 | CoralHue® Humanized Monomeric Azami-Green (phmAG1-S1) | 20 µg |
| AM-V0035 | CoralHue® Humanized Monomeric Azami-Green (phmAG1-MC1) | 20 µg |
| AM-V0036 | CoralHue® Humanized Monomeric Azami-Green (phmAG1-MN1) | 20 µg |
| AM-V0039 | CoralHue® Humanized Monomeric Azami-Green (phmAG1-MCLinker) | 20 µg |
| AM-V0030 | CoralHue® Humanized Monomeric Azami-Green (phmAG1-MNLinker) | 20 µg |
| AM-V0201 | CoralHue® Mitochondria-targeted mAG1 Expression Plasmid (pMT-mAG1) | 20 µg |
| AM-V0202 | CoralHue® ER-targeted mAG1 Expression Plasmid (pER-mAG1) | 20 µg |
| AM-V0203 | CoralHue® Plasma Membrane-targeted mAG1 Expression Plasmid (pPM-mAG1) | 20 µg |
| AM-V0214 | CoralHue® Nucleoplasm-targeted AG Expression Plasmid (pNP-AG) | 20 µg |

| Anti- CoralHue® Azami-Green Antibodies | | | | | | |
|--|--|------------|--------------|--------|--------------|--|
| Code No. | Product | Clone | Isotype | Size | Applications | |
| M102-3 | Anti- CoralHue® Azami Green Monoclonal Antibody | 2F11 | mouse IgG1κ | 100 µg | WB | |
| M102-3S | Anti- CoralHue® Azami Green Monoclonal Antibody (Trial Size) | 2F11 | mouse IgG1κ | 10 µL | WB | |
| M103-3 | Anti- CoralHue® Azami Green Monoclonal Antibody | 3D10 | mouse IgG2ak | 100 µg | IPP | |
| M103-3S | Anti- CoralHue® Azami Green Monoclonal Antibody (Trial Size) | 3D10 | mouse IgG2ak | 10 µL | IPP | |
| PM011 | Anti- CoralHue® Azami Green Polyclonal Antibody | polyclonal | rabbit IgG | 500 µg | WB | |
| PM011S | Anti- CoralHue® Azami Green Polyclonal Antibody (Trial Size) | polyclonal | rabbit IgG | 10 µL | WB | |


Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com

CoralHue[®] Dronpa-Green

CoralHue[®] Dronpa-Green is a monomeric fluorescent protein that has a unique ability to turn on and off its green fluorescence. When subjected to excitation at 400 nm, Dronpa-Green displays a bright green fluorescence which is equally bright as EGFP. When subjected to excitation at 490nm, Dronpa-green's bright green fluorescence is "switched off". Then these bleached proteins can be "switched on" again by exciting the protein at 400nm. This photoconversion can be repeated endlessly, without compromising the brightness of the protein. This unique property of Dronpa-Green is useful for measuring the dynamics of molecular mobility (e.g. diffusion, transport, etc.) of fluorescently labeled molecules in membranes or in living cells.




CoralHue[®] Dronpa-Green is available as several different expression plasmids. Dronpa-Green expression plasmids allow for proteins of interest to be labeled by Dronpa-Green at either their C-terminus or N-terminus. Dronpa-Green is also available as a cDNA plasmid which allows Dronpa-Green to be incorporated into an expression plasmid of choice at several different restriction sites.



CoralHue[®] DG1 expression in HeLa cells

FUN FACT: Dronpa, after "dron" a ninja term for vanishing, and "pa" which stands for photoactivation.

| | Excit./Emiss.Maxima (nm) | Extinction Coefficient(M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|-----|--------------------------|---|----------------------------|----------------|
| DG1 | 503 / 518 | 95,000 (503 nm) | 0.85 | pKa=5.0 |




CoralHue[®] DG1 stained Endoplasmic Reticulum and Actin in HeLa cells.

Photos provided courtesy of Dr. Michael Davidson, National High Magnetic Field Laboratory, The University of Florida.

References

- Fron, E., et al., (2007) J Am Chem Soc. 129, 4870-4871.
- Ando, R., et al., (2007) Biophys J. 92, L97-99.
- Ando, R., et al., (2004) Science 306, 1370-1373.
- Aramaki, S., and Hatta, K., (2006) Dev. Dyn. 235, 2192-2199.
- Dedecker, P., et al., (2006) Biophys. J. 91, 45-47.
- Eisenstein, M., (2005) Nature Methods 2, 8-9.
- Habuchi, S., et al., (2005) PNAS 102, 9511-9516.
- Habuchi, S., et al., (2006) Photochem. Photobiol. Sci. 5, 567-576.
- Kurokawa, K., et al., (2005) Mol. Biol. Cell 16, 4294-4303.
- Sauer, M., (2005) PNAS 102, 9433-9434.
- Stiel, A.C., et al., (2007) Biochem. J. 402, 35-42.
- Wilmann, P.G., et al., (2006) J. Mol. Biol. 364, 213-214.

Vector




| CoralHue ⁺ Dronpa-Green Fluorescent Proteins | | |
|---|---|-------|
| Code No. | Product | Size |
| AM-V0071 | CoralHue ⁺ Dronpa-Green (pDG1-S1) | 20 µg |
| AM-V0072 | CoralHue ⁺ Dronpa-Green (pDG1-MC1) | 20 µg |
| AM-V0073 | CoralHue ⁺ Dronpa-Green (pDG1-MN1) | 20 µg |
| AM-V0131 | CoralHue ⁺ Dronpa-Green3 (pDG3-S1) | 20 µg |

| Anti- CoralHue ⁺ Dronpa-Green Antibodies | | | | | | |
|---|--|-------|-------------|--------|--------------|--|
| Code No. | Product | Clone | Isotype | Size | Applications | |
| M117-3 | Anti- CoralHue ⁺ Dronpa Green Monoclonal Antibody | 4D12 | mouse IgG2a | 100 µL | WB | |
| M118-3 | Anti- CoralHue ⁺ Dronpa Green Monoclonal Antibody | 2F6 | mouse IgG2b | 100 µg | IPP | |

CoralHue[®] Kaede




CoralHue[®] Kaede is a unique fluorescent protein that has the ability to irreversibly convert from bright green fluorescence to bright red fluorescence. This photoconversion can be activated by irradiating the protein with UV or violet light (350-410 nm). Kaede's red fluorescence is comparable in intensity to the protein's green fluorescence and is also stable under usual aerobic conditions, unlike many other photoconverting proteins. The irreversible photoconversion of Kaede provides a simple and powerful technique for regional optical marking.

CoralHue[®] Kaede is available as several different vectors which allow researchers to tag their protein of interest either at the C-terminus or N-terminus. Kaede is also available as a cDNA vector which can be inserted into a plasmid using several different restriction sites.



Kaede expression in HeLa cells
A: Before UV irradiation, B: After UV irradiation *

| | Excit./Emiss.Maxima (nm) | Extinction Coefficient(M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|-------|--------------------------|---|----------------------------|----------------|
| Green | 508 / 518 | 98,800 (508 nm) | 0.88 | pKa=5.6 |
| Red | 572 / 580 | 60,400 (572 nm) | 0.33 | pKa=5.6 |




HeLa cells before and after Kaede photoconversion *

* Photos provided courtesy of Dr. Miyawaki, RIKEN Institute, Japan


References

- Stark, D.A., et al., (2007) Dev. Dyn. 236, 1583-1594.
- Hatta, K., et al., (2006) Nat. Protoc. 1, 960-967.
- Scott, E.K., et al., (2007) Nat. Methods. 4, 323-326.
- Davison, J.M., et al., (2007) Dev. Biol. 304, 811-824.
- Pisharath, H., et al., (2007) Mech. Dev. 124, 218-229.
- Ando, R., et al., (2002) PNAS 99, 12651-12656.
- Arimura, S., et al., (2004) PNAS 101, 7805-7808.
- Dittrich, P.S., et al., (2005) Biophysical J. 89, 3446-3455.
- Hosoi, H., et al., (2006) J. Phys. Chem. B. 110, 22853-22860.
- Kimura, Y., et al., (2005) J. Neurosci. 26, 5684-5697.
- Kurokawa, K., et al., (2005) Mol. Biol. Cell 16, 4294-4303.
- Lippincott-Schwartz, J., et al., (2003) Imaging in Cell Biology, 57-514.
- Melton, L., et al., (2005) NATURE 437, 775-779.
- Miyawaki, A., (2002) Cell Structure and Function 27, 343-347.
- Miyawaki, A., (2004) Nature Biotechnology 22, 1374-1376.
- Mizuno, H., et al., (2003) Molecular Cell 12, 1051-1058.
- Mutoh, T., et al., (2006) Exp. Neurol. 200, 430-437.
- Raab-Graham, K.F., et al., (2006) Science 314, 144-148.
- Sato, T., et al., (2006) Genesis 44, 136-142.
- Schafer, S.P., et al., (2006) Microsc. Res. Tech. 69, 210-219.
- Tsutsui, H., et al., (2005) EMBO reports 6, 1-6.
- Wiedenmann, J., et al., (2004) PNAS 101, 15905-15910.



Neurons transfected with Kaede (before and after photoconversion) *

Vector



| CoralHue [®] Kaede Fluorescent Proteins | | |
|--|--|-------|
| Code No. | Product | Size |
| AM-V0011 | CoralHue [®] Kaede (pKaede-S1) | 20 µg |
| AM-V0012 | CoralHue [®] Kaede (pKaede-MC1) | 20 µg |
| AM-V0013 | CoralHue [®] Kaede (pKaede-MN1) | 20 µg |

| Anti-CoralHue [®] Kaede Antibodies | | | | | |
|---|--|------------|-------------|-------|--------------|
| Code No. | Product | Clone | Isotype | Size | Applications |
| M106-3 | Anti-CoralHue [®] Kaede Monoclonal Antibody | 2F4 | mouse IgG1k | 100µg | IPP |
| M106-3S | Anti-CoralHue [®] Kaede Monoclonal Antibody(Trial Size) | 2F4 | mouse IgG1k | 10µL | IPP |
| M125-3 | Anti-CoralHue [®] Kaede Monoclonal Antibody | 3B1 | mouse IgG1 | 100µL | WB |
| PM012 | Anti-CoralHue [®] Kaede Polyclonal Antibody | polyclonal | rabbit IgG | 500µL | IPP |
| PM012S | Anti-CoralHue [®] Kaede Polyclonal Antibody(Trial Size) | polyclonal | rabbit IgG | 10µL | WB |

Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com

CoralHue[®] Keima-Red

MBL International's new, exclusive fluorescent proteins, CoralHue[®] Monomeric and Dimeric Keima-Red, combine bright red fluorescence with the largest commercially available Stokes shift (ex. 440 nm, em. 620 nm), making Keima-Red a superb reporter protein for multicolor fluorescence analyses. Keima-Red is particularly useful when performing dual-color fluorescence cross-correlation spectroscopy (FCCS) because it can be paired with a fluorescent protein with a similar excitation wavelength that has a small Stokes shift. This pairing would allow for simultaneous excitation of the two proteins without interference between the two emissions.






CoralHue[®] Keima-Red expression in HeLa cells


FUN FACT: Keima means "knight" of Japanese chess (shogi). The protein named so for its "jumping" Stokes shift.

The CoralHue[®] Keima-Red protein cloning plasmids allow for insertion of cDNA sequences to create protein fusion products between the protein of interest and Keima. One insertion locations is available allowing the target protein to be tagged by Keima at N-terminus. CoralHue[®] Keima-Red fluorescent protein cloning plasmids create protein fusion products that are useful for tracking protein localization within cells as well as monitoring gene expression. Keima-Red is also now available as target specific constructs which will allow for Keima-Red protein fusion products to be directed to either the mitochondria or the plasma membrane.

| | Excit./Emiss.Maxima (nm) | Extinction Coefficient (M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|------------|--------------------------|--|----------------------------|----------------|
| dKeima-Red | 440 / 616 | 24,600 (at 440 nm) | 0.31 | pKa=6.5 |
| mKeima-Red | 440 / 620 | 14,400 (at 440 nm) | 0.24 | pKa=6.5 |




CoralHue[®] MT-Keima-Red expression in HeLa cells



CoralHue[®] PM-Keima-Red expression in HeLa cells

References

Kogure, T., et al., (2006) Nat. Biotechnol. 24, 577-581.




CoralHue[®] mKeima-Red
Zyxin

Golgi

Keratin

Photos provided courtesy of Dr. Michael Davidson, National High Magnetic Field Laboratory, The University of Florida.

Vector



| CoralHue [®] Keima-Red Fluorescent Proteins | | |
|--|---|-------|
| Code No. | Product | Size |
| AM-V0091 | CoralHue [®] Monomeric Keima-Red (pmKeima-Red-S1) | 20 µg |
| AM-V0093 | CoralHue [®] Monomeric Keima-Red (pmKeima-Red-MN1) | 20 µg |
| AM-V0101 | CoralHue [®] Dimeric Keima-Red (pdKeima-Red-S1) | 20 µg |
| AM-V0121 | CoralHue [®] Dimeric Keima570 (pdKeima570-S1) | 20 µg |
| AM-V0251 | CoralHue [®] Mitochondria-targeted mKeima-Red Expression Plasmid (pMT-mKeima-Red) | 20 µg |
| AM-V0253 | CoralHue [®] Plasma Membrane-targeted mKeima-Red Expression Plasmid (pPM-mKeima-Red) | 20 µg |

| Anti- CoralHue [®] Keima-Red Antibodies | | | | | | |
|--|---|-------|-------------|--------|--------------|--|
| Code No. | Product | Clone | Isotype | Size | Applications | |
| M126-3 | Anti- CoralHue [®] Keima-Red Monoclonal Antibody | 2F7 | mouse IgG2a | 100 µg | WB | |
| M127-3 | Anti- CoralHue [®] Keima-Red Monoclonal Antibody | 3C9 | mouse IgG1 | 100 µg | IPP | |


Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com

CoralHue[®] Kikume Green-Red




CoralHue[®] Kikume Green-Red (KikGR) is a photoconverting protein which has the capability to irreversibly convert from green to red fluorescence when subjected to UV or violet light. The excitation lights used to elicit red or green fluorescence do not induce the photoconversion. This provides a simple and powerful technique for regional optical marking. KikGR's red and green fluorescence can be activated in vivo and is several-fold brighter than Kaede's green and red fluorescence.

CoralHue[®] KikGR is available as several different plasmids allowing for several different insertion sites, including the N-terminus and the C-terminus of a protein of interest. CoralHue[®] KikGR is now also available as a humanized plasmid which can be expressed in mammalian cells. CoralHue[®] KikGR's ability to photoconvert from green to red combined with the many forms of the protein available make KikGR the perfect fluorescent protein for regional optical marking.



CoralHue[®] hKikGR1 expression in HeLa cells. The fluorescence of hKikGR1 irradiated with UV is red. Untreated hKikGR emits green fluorescence.

| | Excit./Emiss.Maxima (nm) | Extinction Coefficient(M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH Sensitivity |
|-------|--------------------------|---|----------------------------|----------------|
| Green | 507 / 517 | 53,700 (507 nm) | 0.70 | pKa=7.8 |
| Red | 583 / 593 | 35,100 (583 nm) | 0.65 | pKa=5.5 |




CoralHue[®] KikGR1 Mitochondria Endoplasmic reticulum

The fluorescence of KikGR1 irradiated with UV is red. Untreated KikGR emits green fluorescence. Photos provided courtesy of Dr. Michael Davidson, National High Magnetic Field Laboratory, The University of Florida.

Reference
 Stark, D.A., et al., (2007) Dev Dyn. 236, 1583-1594.
 Tsutsui, H., et al., (2005) EMBO reports 6, 1-6.

Vector



| CoralHue [®] Kikume Green-Red Fluorescent Proteins | | |
|---|--|-------|
| Code No. | Product | Size |
| AM-V0081 | CoralHue [®] Kikume Green-Red (pKikGR1-S1) | 20 µg |
| AM-V0082 | CoralHue [®] Kikume Green-Red (pKikGR1-MC1) | 20 µg |
| AM-V0083 | CoralHue [®] Kikume Green-Red (pKikGR1-MN1) | 20 µg |
| AM-V0084 | CoralHue [®] Humanized Kikume Green-Red (phKikGR1-S1) | 20 µg |
| AM-V0085 | CoralHue [®] Humanized Kikume Green-Red (phKikGR1-MC1) | 20 µg |
| AM-V0086 | CoralHue [®] Humanized Kikume Green-Red (phKikGR1-MN1) | 20 µg |
| AM-V0089 | CoralHue [®] Humanized Kikume Green-Red (phKikGR1-MCLinker) | 20 µg |
| AM-V0080 | CoralHue [®] Humanized Kikume Green-Red (phKikGR1-MNLinker) | 20 µg |

| Anti-CoralHue [®] Kikume Green-Red Antibodies | | | | | |
|--|---|-------|-------------|--------|--------------|
| Code No. | Product | Clone | Isotype | Size | Applications |
| M128-3 | Anti-CoralHue [®] Kikume Green-Red/KikGR Monoclonal Antibody | 5B3 | mouse IgG2b | 100 µg | WB |
| M129-3 | Anti-CoralHue [®] Kikume Green-Red/KikGR Monoclonal Antibody | 2D3 | mouse IgG2b | 100 µg | IPP |


Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com

CoralHue[®] Kusabira-Orange




CoralHue[®] Kusabira-Orange is a fluorescent protein that was derived from the stony coral whose Japanese name is “Kusabira-ishi”. Kusabira-Orange absorbs light maximally at 548 nm and emits orange light at 561 nm. Wild-Type Kusabira-Orange rapidly matures to form a fluorescent dimeric complex which can be used to mark cells or to report gene expression without problems stemming from protein aggregation. CoralHue[®] Kusabira-Orange has also been engineered as monomeric and humanized forms.

CoralHue[®] Kusabira-Orange is available as plasmids which are fusion to the C-terminus or the N-terminus of your proteins of interest. Kusabira-Orange is available as several targeted expression plasmids that are specific to the endoplasmic reticulum, the nucleoplasm, and the plasma membrane, and the mitochondria.




CoralHue[®] KO1 expressed in olfactory neurons in fish. Photo provided courtesy of Dr. Yoshihara, RIKEN Institute, Japan.


| | Excit. /Emiss.Maxima (nm) | Extinction Coefficient(M ¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|------|---------------------------|--|----------------------------|-----------------------|
| mKO1 | 548 / 559 | 51,600 (548 nm) | 0.6 | pK _a = 5.0 |



Kusabira Orange Plasma Membrane



Kusabira Orange Vimentin *



Kusabira Orange Mitochondria *




Kusabira Orange Nucleoplasm

* Photo provided courtesy of Dr. Michael W. Davidson, The National High Magnetic Field Laboratory, The Florida State University

References

Niwa, H., et. al., (2005) Cell 123, 917-929.
 Shaner, N.C., et. al., (2005) Nat. Methods. 2, 905-909. Review.
 Karasawa, S., et. al., (2004) Biochem. J. 381, 307-312.
 Ishida, A., et. al., (2005) Res. Bull. Aichi Agric. Res. Ctr. 37, 141-146.

Vector



| CoralHue [®] Kusabira-Orange Fluorescent Proteins | | |
|--|---|--------|
| Code No. | Product | Size |
| AM-V0041 | CoralHue [®] Kusabira-Orange (pKO1-S1) | 2.0 µg |
| AM-V0051 | CoralHue [®] Monomeric Kusabira-Orange (pmKO1-S1) | 2.0 µg |
| AM-V0052 | CoralHue [®] Monomeric Kusabira-Orange (pmKO1-MC1) | 2.0 µg |
| AM-V0053 | CoralHue [®] Monomeric Kusabira-Orange (pmKO1-MN1) | 2.0 µg |
| AM-V0044 | CoralHue [®] Humanized Kusabira-Orange (phKO1-S1) | 2.0 µg |
| AM-V0045 | CoralHue [®] Humanized Kusabira-Orange (phKO1-MC1) | 2.0 µg |
| AM-V0046 | CoralHue [®] Humanized Kusabira-Orange (phKO1-MN1) | 2.0 µg |
| AM-V0054 | CoralHue [®] Humanized Monomeric Kusabira-Orange (phmKO1-S1) | 2.0 µg |
| AM-V0055 | CoralHue [®] Humanized Monomeric Kusabira-Orange (phmKO1-MC1) | 2.0 µg |
| AM-V0056 | CoralHue [®] Humanized Monomeric Kusabira-Orange (phmKO1-MN1) | 2.0 µg |
| AM-V0059 | CoralHue [®] Humanized Monomeric Kusabira-Orange (phmKO1-MCLinker) | 2.0 µg |
| AM-V0050 | CoralHue [®] Humanized Monomeric Kusabira-Orange (phmKO1-MNLinker) | 2.0 µg |
| AM-V0221 | CoralHue [®] Mitochondria-targeted mKO1 Expression Plasmid (pMT-mKO1) | 2.0 µg |
| AM-V0222 | CoralHue [®] ER-targeted mKO1 Expression Plasmid (pER-mKO1) | 2.0 µg |
| AM-V0223 | CoralHue [®] Plasma Membrane-targeted mKO1 Expression Plasmid (pPM-mKO1) | 2.0 µg |
| AM-V0234 | CoralHue [®] Nucleoplasm-targeted KO Expression Plasmid (pNP-KO) | 2.0 µg |


| Anti- CoralHue [®] Kusabira-Orange Antibodies | | | | | |
|--|--|-------|-------------|--------|--------------|
| Code No. | Product | Clone | Isotype | Size | Applications |
| M104-3 | Anti- CoralHue [®] Kusabira Orange Monoclonal Antibody | 1H7 | mouse IgG1κ | 100 µg | WB |
| M104-3S | Anti- CoralHue [®] Kusabira Orange Monoclonal Antibody (Trial Size) | 1H7 | mouse IgG1κ | 10 µL | WB |
| M105-3 | Anti- CoralHue [®] Kusabira Orange Monoclonal Antibody | 2G9 | mouse IgG1κ | 100 µg | IPP |
| M105-3S | Anti- CoralHue [®] Kusabira Orange Monoclonal Antibody (Trial Size) | 2G9 | mouse IgG1κ | 10 µL | IPP |

Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com



CoralHue[®] Midoriishi-Cyan

CoralHue[®] Midoriishi-Cyan fluorescent protein is derived from a stony coral who's Japanese name is "Midori-ishi". Midoriishi-Cyan absorbs light maximally at 472 nm and emits cyan light at 495 nm. Wild-type Midoriishi-Cyan rapidly matures to form a fluorescent dimeric complex which can be used to mark individual cells or to report gene expression without problems stemming from protein aggregation.




CoralHue[®] MiCy1 expression in HeLa cell

| | Excit. /Emiss.Maxima (nm) | Extinction Coefficient(M ⁻¹ cm ⁻¹) | Fluorescence Quantum Yield | pH sensitivity |
|-------|---------------------------|--|----------------------------|----------------|
| MiCy1 | 472 /495 | 27,250 (472 nm) | 0.9 | pK a = 6. 6 |



Vector



References
 Shaner, N.C., et al., (2005) Nat. Methods. 2, 905-909. Review.
 Karasawa, S., et al., (2004) Biochem. J. 381, 307-312.

| CoralHue [®] Midoriishi-Cyan Fluorescent Proteins | | |
|--|---|-------|
| Code No. | Product | Size |
| AM-V0061 | CoralHue [®] Midoriishi-Cyan (pMiCy1-S1) | 20 µg |

| Anti- CoralHue [®] Midoriishi-Cyan Antibodies | | | | | |
|--|---|-------|-------------|--------|--------------|
| Code No. | Product | Clone | Isotype | Size | Applications |
| M116-3 | Anti- CoralHue [®] Midoriishi-Cyan Monoclonal Antibody | 2C1 | mouse IgG2b | 100 µg | IPP |
| M130-3 | Anti- CoralHue [®] Midoriishi-Cyan Monoclonal Antibody | 5B7 | mouse IgG1 | 100 µg | WB |

Application; WB: Western blotting, IPP: Immunoprecipitation

For more information, go to www.mblintl.com


Targeted Plasmids

Azami-Green


Kusabira-Orange

Keima-Red

Plasma membrane Targeting



CoralHue[®] pPM-mAG1




CoralHue[®] pPM-mKO1




CoralHue[®] pPM-mKeima-Red


Mitochondria Targeting



CoralHue[®] pMT-mAG1




CoralHue[®] pMT-mKO1




CoralHue[®] pMT-mKeima-Red

Endoplasmic reticulum Targeting




CoralHue[®] pER-mAG1



CoralHue[®] pER-mKO1


Nucleoplasm Targeting



CoralHue[®] pNP-AG



CoralHue[®] pNP-KO



| Code No. | Product | Size |
|----------|---|-------|
| AM-V0203 | CoralHue [®] Plasma Membrane-targeted mAG1 Expression Plasmid (pPM-mAG1) | 20 µg |
| AM-V0223 | CoralHue [®] Plasma Membrane-targeted mKO1 Expression Plasmid (pPM-mKO1) | 20 µg |
| AM-V0253 | CoralHue [®] Plasma Membrane-targeted mKeima-Red Expression Plasmid (pPM-mKeima-Red) | 20 µg |
| AM-V0201 | CoralHue [®] Mitochondria-targeted mAG1 Expression Plasmid (pMT-mAG1) | 20 µg |
| AM-V0221 | CoralHue [®] Mitochondria-targeted mKO1 Expression Plasmid (pMT-mKO1) | 20 µg |
| AM-V0251 | CoralHue [®] Mitochondria-targeted mKeima-Red Expression Plasmid (pMT-mKeima-Red) | 20 µg |
| AM-V0202 | CoralHue [®] ER-targeted mAG1 Expression Plasmid (pER-mAG1) | 20 µg |
| AM-V0222 | CoralHue [®] ER-targeted mKO1 Expression Plasmid (pER-mKO1) | 20 µg |
| AM-V0214 | CoralHue [®] Nucleoplasm-targeted AG Expression Plasmid (pNP-AG) | 20 µg |
| AM-V0234 | CoralHue [®] Nucleoplasm-targeted KO Expression Plasmid (pNP-KO) | 20 µg |

CoralHue[®] fluorescent proteins used in these products were co-developed with the Laboratory for Cell Function and Dynamics, the Advanced Technology Development Center, the Brain Science Institute, and the Institute of Physical and Chemical Research (RIKEN) (lab head Dr. Atsushi Miyawaki).

The use of CoralHue[®] fluorescent proteins require a license. Amalgaam grants non-profit research organizations an international, royalty-free, non-exclusive, limited license to use this product for non-commercial research use only. This license excludes the right to sell or transfer this product, its components, or modifications of this product to third parties. Any other uses require a license. The use of this product by for profit organizations, for either commercial or non-commercial use, requires a license. For additional licensing information, please contact MBL at 781-939-6964 or sales@mblintl.com. The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242."

Vector sequences are available on our website under Technical Resources.

MBL International Corporation is a leading life science company focused on providing high quality products and solutions for life science research and clinical diagnostics.

MBL International Corporation develops, manufactures and markets a wide range of ELISA and gold standard IFA test kits to aid in the diagnosis of autoimmune and infectious disease as well as an extensive portfolio of monoclonal and polyclonal antibodies, fluorescent and recombinant proteins, and ELISA kits for research related to cancer, immunology, neuroscience, allergy, autophagy, and apoptosis.

Our products are used widely in academic research institutions, pharmaceutical and biotechnology companies, government agencies, as well hospital and reference laboratories. By providing a consultative approach, superior technical and customer support, and convenient purchasing options, MBL International Corporation is your partner in solving your scientific and clinical challenges.

MBL
International Corporation
4 H Constitution Way
Woburn, MA 01801, USA

Tel: (781) 939-6964
(800) 200-5459
Fax: (781) 939-6963

URL: <http://www.mblintl.com> email: info@mblintl.com

MBL
Medical & Biological Laboratories Co., LTD.
KDX Nagoya Sakae Bldg. 10F
4-5-3 Sakae, Naka-ku, Nagoya, Aichi
460-0008, JAPAN
Tel: 81-52-971-2081 Fax: 81-52-971-2337
URL: <https://ruo.mbl.co.jp>