

# ExoCap™ Ultracentrifugation/Storage Booster

Code No.

Quantity

MEX-USB

50 mL

## Product description

Exosomes are extracellular vesicles secreted by most cell types and contain various marker proteins and RNAs, such as microRNA and fragmented mRNA. **ExoCap™ Ultracentrifugation/Storage Booster (ExoCap™ USB)** is an additive for exosome-conditioned media specifically designed to 1) improve the recovery rate of ultracentrifugation and 2) stabilize and extend 4°C exosome storage. The advantages of this reagent are as follows:

- Improving exosome recovery rate from cell culture supernatant by ultracentrifugation
- Improving stability at 4°C storage for exosomes purified from cell culture supernatant
- Acceptable at -80°C storage for exosome
- No complicated troublesome work

## Product contents

Component	Amount	Storage
ExoCap™ Ultracentrifugation/Storage Booster	50 mL	2-8°C

**Note:** This reagent contains no preservative; handle under aseptic conditions.

## Storage and Stability

Stable for 2.5 years from the date of receipt when stored at 2-8°C. Do not freeze.

**Note: If precipitation occurs during storage, equilibrate the reagent to room temperature and mix well to re-dissolve the precipitation.**


## Materials required but Not Provided

1. Refrigerator (2°C to 8°C)
2. Ultra low temperature freezer (-80°C)
3. Ultracentrifuge capable of 100,000 × g  
(e.g., Optima XE-90, Rotor: SW32Ti and SW40Ti (Beckman Coulter))
4. Ultracentrifuge tubes (e.g., Beckman Coulter; code no. 344058 and 344060)
5. Centrifuge tubes [Recommendation; use low-adhesion tube]
6. Pipette tips [Recommendation; use low-adhesion pipette tip]
7. Pipettes (Micropipettes)
8. 0.22 µm filter
9. MilliQ water

## ExoCap™ Ultracentrifugation/Storage Booster Procedure

### A. Ultracentrifugation for exosome isolation

#### 1. Procedure summary



#### 2. Buffer preparation

**Recovery buffer** (for the 2<sup>nd</sup> ultracentrifugation and the final re-dispersion of exosomes)

Mix 1 part of the **ExoCap™ USB** with 9 parts of MilliQ water before use.

Adjust the Recovery buffer preparation according to the size of ultracentrifugation instruments and apparatuses.

#### 3. Protocols

##### ◆ Pre-clearance

1. Collect the serum-free exosome-conditioned cultured medium from the cell culture dish into the centrifugation tube.
2. Centrifuge the tube at  $300 \times g$  at  $4^\circ\text{C}$  for 10 minutes.
3. Transfer the supernatant to a new centrifugation tube and discard cell pellet.
4. Centrifuge the tube at  $2,000 \times g$  at  $4^\circ\text{C}$  for 10 minutes.
5. Transfer the supernatant to a new tube and discard dead cell pellet.
6. Filter the final supernatant with  $0.22 \mu\text{m}$  filter.

##### ◆ Ultracentrifugation and storage

The required volume of ExoCap™ USB is described in the table below.

Sample	ExoCap™ USB	Recovery Buffer		Total use of ExoCap™ USB /sample
	1 <sup>st</sup> UC	2 <sup>nd</sup> UC <sup>†</sup>	Storage <sup>†</sup>	
10 mL	1.1 mL	10 mL (1.1 mL)	0.1 mL (0.01 mL)	2.21 mL
30 mL	3.3 mL	30 mL (3.3 mL)	0.3 mL (0.03 mL)	6.63 mL
50 mL	5.5 mL	50 mL (5.5 mL)	0.5 mL (0.55 mL)	11.65 mL

<sup>†</sup>Numbers in parentheses shows the amount of ExoCap™ USB.

**Note:** The usage depends on the ultracentrifugation scale to use. When smaller size tube is used for the 2<sup>nd</sup> ultracentrifugation (2<sup>nd</sup> UC), the total amount of ExoCap™ USB can be decreased.

1. Add 1 part of the **ExoCap™ USB** to 9 parts of the supernatant (prepared by pre-clearance step) and mix well.  
(Example) Supernatant 30 mL + ExoCap™ USB 3.3 mL

2. Transfer the sample to the ultracentrifugation tube (*e.g.* Beckman Coulter; code no. 344058).
3. 1<sup>st</sup> Ultracentrifugation.  
(Example) Apparatus: Optima XE-90, Rotor: SW32Ti (Beckman Coulter), 100,000 × *g* at 4°C for 70 minutes.  
**Note:** During the 1<sup>st</sup> ultracentrifugation, prepare the **Recovery buffer**.  
(Example) ExoCap™ USB 1.2 mL + MilliQ water 10.8 mL
4. Discard the supernatant.
5. Suspend the pellet by Recovery buffer and transfer it to the new ultracentrifugation tube (*e.g.* Beckman Coulter; code no. 344060).  
(Example) The pellet is suspended by 11 mL of **Recovery buffer**. Transfer it to the new ultracentrifugation tube.
6. 2<sup>nd</sup> ultracentrifugation.  
(Example) Apparatus: Optima XE-90, Rotor: SW40Ti (Beckman Coulter), 100,000 × *g* at 4°C for 70 minutes.
7. Discard the supernatant.
8. Suspend the pellet with **Recovery buffer**.  
(Example) The pellet is resuspended with 0.3 mL **Recovery buffer** (100 times concentrated).
9. Storage at 4°C or -80°C in portions.

## **B. Storage of isolated exosomes**

### **Protocols**

1. Purify exosome from cell culture supernatant by proper method such as ultracentrifugation, size exclusion chromatography, or Beads.
2. Disperse the exosomes in proper buffer (*e.g.* PBS).
3. Add 1 part of the ExoCap™ USB to 9 parts of the sample and mix well.
4. Storage at 4°C or -80°C in portions.

**Note:** Even if ExoCap™ USB is used for storage, repeated freeze-thaw or storage at 4°C will occasionally cause a deterioration of the collected exosomes. Therefore, for long-term storage, preparation of small aliquots for -80°C storage is recommended. The exosome solution with the additive at 4°C should be used within a few days.


### **Attention**

- The efficiency of ExoCap™ USB depends on the sample. Preliminary tests are strongly recommended.
- According to the purpose of use, prepared exosome solution should be sterilized before use.

## Examples of ExoCap™ USB


### Exosome Purification by Ultracentrifugation

#### 1) Evaluation of the ExoCap™ USB effect for ultracentrifugation by Western blotting (Non-reduced)



HT29 cells were cultured up to 60-80%. Subsequently, wash the cells and change the new medium without FBS. 3 days later, recover the supernatant. The cell density was reached from 90% to confluent. After the pre-clearance by centrifugation at  $300-3,000 \times g$ , this supernatant was filtrated by  $0.22 \mu\text{m}$  filter and then used for the experiment. Exosome was isolated from the supernatant by ultracentrifugation with/without ExoCap™ USB. The intensity of the bands of exosome markers with the ExoCap™ USB was denser than those of without ExoCap™ USB. Exosome recovery was boosted by using the ExoCap™ USB.

#### 2) Evaluation of the ExoCap™ USB effect for ultracentrifugation by CLEIA



Exosome recovery rate was improved as a result of the ExoCap USB addition. For beads CLEIA (Chemiluminescent Enzyme Immuno Assay) ExoCap™ Streptavidin Kits, biotinylated CD81 antibody and alkaline phosphatase (ALP) labeled CD9 antibody were used.

## Manufacturer

**MBL** MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL <http://ruo.mbl.co.jp>

e-mail [support@mbi.co.jp](mailto:support@mbi.co.jp), TEL 052-238-1904